

How-to: Use Syringe Pump for line flushes



First:
Turn on (counter clock-wise) the compressed air in the O-18 fume hood.



The regulator should be set to **2 bar** when the air is flowing.



When the “off” lever is pointed away from the blue line attached to the pneumatic cylinder, the air flow is blocked and the plunger holder should move freely up and down.



When the “off” lever is pointed towards the empty connection (pressure release) the air flow is no longer restricted and air flows into the pneumatic cylinder which causes the plunger holder to constrict towards the cylinder body

How-to: Use Syringe Pump for line flushes



Attach syringe(s) to final product lines



Place syringes in their designated slots prior to air open. The flanges should be flush against the flange holder, and no gaps should be present (see red arrow on left syringe).



"Off" lever is restricting the air flow so syringes can be placed in correct position.

Note: When changing lever position there should be an audible hiss when alleviating pressure from the lines. This is normal as the compressed air from inside the syringe needs to be exhausted.



"Off" lever is in middle position and the cylinder constricts. Make sure full contents of syringe are dispensed.

CAUTION: Hold plunger upright as much as possible as not to damage the lines once the syringe pump is installed.